

Forever Medical®

MRI ANESTHESIA

Anesthesia for MRI





YJ-A800



YJ-A801



YJ-A802



YJ-A803




YJ-A804



YJ-A805

Model	YJ-A800	YJ-A801	YJ-A802	YJ-A803	YJ-A804	YJ-A805
Part of Anesthesia Machine						
Main body	High-strength engineering plastic rack, light, beautiful and corrosion resistant					
Scope of application	adult, child(4 years above)					
Gas source	O2:0.28~0.6 MPa		O2, N2O: 0.28~0.6MPa			O2, N2O, Air: 0.28~0.6M Pa
Flow meter	O2:	O2:	O2,N2O:	O2,N2O:	O2,N2O:	O2,N2O:
	n1.0-10L/min	0.051.0L/min 1.1-10L/min	0.0510L/min	0.05-1.0L/min	0.05-1.0L/min	0.05-1.0L/min
			0.05-10L/min	1.1-10L/min	1.1-10L/min	1.1-10L/min
						Air:0.05-10L/min
O2,N2O linkage and N2O Stopper	/	/	When using nitrous oxide, oxygen concentration>25%. When the oxygen pressure<0.2MPa, the flow of nitrous oxide would be cut off.			
Vaporizer	one	One, high-precision evaporator	One, high-precision evaporator	One, high-precision evaporator	Two, high-precision evaporator	Two,high-precision evaporator
Part of Ventilator						
Display mode	/	5.7"	5.7"	10.4" TFT	5.7"	10.4" TFT
mode	/	IPPV, SIPPV, SIMV, MANUAL	IPPV,SIPPV, SIMV, MANUAL	IPPV, SIPPV, SIMV,MANUAL, VCV, PCV	IPPV,SIPPV, SIMV, MANUAL	IPPV, SIPPV, SIMV,MANUAL, VCV, PCV
function	/	PEEP, SIGH, IRV	PEEP, SIGH, IRV	PEEP,SIGH, IRV, IP, standby	PEEP, SIGH, IRV	PEEP,SIGH, IRV, IP, standby
TV	/	50~1500ml	50~1500ml	20~1600 ml	50~1500ml	20~1600 ml
Rate	/	2~99 bpm	2~99 bpm	2~99 bpm	2~99 bpm	2~99 bpm
SIMV rate	/	2~20 bpm	2~20 bpm	2~20 bpm	2~20 bpm	2~20 bpm
I:E	/	2:1~1:8	2:1~1:8	4:1~1:8	2:1~1:8	4:1~1:8
Inspiratory trigger pressure	/	10~20cmH2O	-10~20cmH2O	-10~20cmH2O	-10~20cmH2O	-10~20cmH2O
PEEP	/	2~20cmH2O	2~20cmH2O	2~20cmH2O	2~20cmH2O	2~20cmH2O
Pressure limitation	/	10~60cmH2O		5~60cmH2O	10~60cmH2O	5~60cmH2O
SIGH	/	1.5 times the inspiratory time/100bpm		SIGH 1.5 times the inspiratory time (60-120 adjustable)	1.5 times the inspiratory time/100bpm	IGH 1.5 times the inspiratory time (60-120 adjustable)
Inspiratory platform	/	/	/	0~50%	/	0~50%
Pressure control	/	/	/	3~50cmH2O	/	3~50cmH2O
Pressure slope	/	/	/	1~10 gear	/	1~10 gear
Power	/	AC 220V 50Hz (UPS, Storage battery)				
Parameters for ventilation monitoring	/	801, 802, 804: tidal volume, minute ventilation volume, IPPV rate, SIMV rate, total respiratory rate, I/E, peak pressure of airway, pressure - time waveform, flow - time waveform, PEEP, inspiratory trigger pressure				
		B3, B5: tidal volume, minute ventilation volume, IPPV rate, SIMV rate, total respiratory rate, I/E, peak pressure of airway, pressure - time waveform, flow rate - time waveform, PEEP, inspiratory trigger pressure, inhale tidal volume, inspiratory platform, average pressure, autonomous inspiratory rate, oxygen concentration, battery capacity, flow-volume loop, pressure-volume loop				
Mechanical Arms	/	/	/	YES	YES	YES
End-tidal CO2 monitoring	Optional					



About us

Henan Forever Medical CO.,LTD is located in the capital of Henan International City- Zhengzhou, enterprise since its establishment in the market practice, trial and error learning, and accumulated a great deal of entrepreneurial experience, forming a collection development, production, processing and marketing in one of entity.

Forever Medical is a company specialized in materials and related products for the chemical composition, physical performance analysis, testing instruments and equipment suppliers and service providers. The products series are: Lab Equipment, ICU Equipment, Medical Imaging, Cold Chain Medical deep freezer, Medical consumables. Mainly products are : Anesthesia System Machine, Patient Monitor, Ventilator, Ultrasound System Machine, Hematology and Biochemistry Analyzer, X-ray System, Dental Sterilization Machine, Surgical Pliers, and Hospital Patient Bed.

Forever Medical®

Henan Forever Medical CO.,LTD

Sales Office: Room 806, Floor 8, Wanda Center, Jinshui District, Zhengzhou, 450000, Henan, China

Tel: +86-371-61315851

Fax: +86-371-61315853

Factory: Industrial area, Changyuan Town, Xinxiang City, Henan Province, 450000, China

Phone: +86-156-37181323 +86-156-37181685

E-mail: forermed@yjcompany.com

Website: www.yongjiumed.com

