

Good Sales Patient Monitor Vital Signs Monitor with 12.1 LED Screen



Introduction of the patient monitor YJ- F6:

The YJ-F6 patient monitor is designed to match the pace and unique needs of adult, pediatric, and neonatal intensive care; anesthesia and peri-operative care; OR and cardiac care environments. Pacemaker detection and Arrhythmia analysis. An optional USB Mouse/Keyboard control eases navigation and data entry. Humane fashion design and handle of humanized integration design, easy to carry. Designed without fan, characterized with no noise and lower consumption.

Key points: We produce Patient Monitor price is cheapest and the patient monitor used in hospital with to monitor vital sign of patient, like ECG, PR, SpO2, Resp., emp. NIBP and so on. It can be used in normal patient room, ICU, NICU.

Good Sales Patient Monitor Vital Signs Monitor with 12.1 LED Screen

Specifications of patient monitor YJ- F6:

Standard Config:	12.1 inch LED, 3/5 Lead ECG, NIBP, SpO2, Pulse Rate, Temperature, Respiration
Optional Config:	2-Temperature, 1/2 IBP, Mainstream EtCO2, Nellcor SpO2, Masimo SpO2
Optional Accessories:	Touch Screen, Printer, Central Monitoring Station, Neo/Ped Accessories, Wall-mounting, Trolley

Production feature and application of patient monitor YJ- F6:

The patient monitor YJ- F6 applied to adult, pediatric, and neonatal intensive care; anesthesia and peri-operative care; or and cardiac care environments. Pacemaker detection and Arrhythmia analysis.

Newly advanced A9 main board with Linux OS, calculating speed is 4 times faster than traditional products. Support storage of 2160 hours trend table and graph review, 2 hours waveform review,

2000 groups NIBP review and 2000 alarm events review.

Production Detail of patient monitor YJ- F6:

Body	12.1 inch high-brightness TFT LED
	Support display 9~12 waveforms
	Support 7 channel ECG waveform display
	Optional touch screen
Printer	Built-in high-speed 50mm thermal printer (Brand: SEIKO, Japan)
Central System	Support wired/wireless connection to the central monitoring station.
Alarm	Three-level acousto-optic alarm
	Sensor-off alarm
	Paper out alarm
	Support alarm review
	Support alarm pause
Operating System	Support operation with mouse and computer keyboard (option)
	Multi- display mode to meet different clinical requirement, including standard, big font, Oxy CRG, trend graph, NIBP review and full leads ECG interface.
	NIBP self-test mode: including over pressure test, static pressure test and air leakage test.
	Generate ID automatically when register a new patient
	Support medical history search by patient ID, name and mobile number
	18 types of Arrhythmia analysis and real-time S-T segment analysis and pacemaker detection
	Drug calculation and tit-ratation table
	Support multi- language display
	Support online software upgrading by net/USB
Interface	Transducer socket
	ECG
	SPO2
	NIBP
	IBP 1 / IBP 2
	TEMP 1 / TEMP 2 / Mainstream EtCO2
LAN / Central monitor system interface	The LAN port is for connecting the monitor to a central monitor system network.
Input device interface	Knob (Standard)
	USB interface (optional)
	Touch screen (optional)

This optional interface provides an USB ports to enable the monitor to be connected to

off-the-shelf input devices:

Mouse: any specified trackball or USB mouse may be used for navigation and data entry.

Computer keyboard: an USB computer keyboard can be used for data entry instead of the on-screen pop-up keyboard.

[Remark: user has to restart the machine after plug and unplug the mouse / computer keyboard.]

· Knob (Standard)

Service Features of the patient monitor YJ- F6:

The Support Tool helps technical personnel to carry out configuration, upgrades and troubleshooting on an individual monitor.

back up the monitor settings.

The Service Mode is password-protected and ensures that only trained staff can access service tests and tasks. It includes Skin Type, Brightness Setting, Clear Data and NIBP Calibration

Production qualification of the patient monitor YJ- F6:

CE, ISO certification

Delivery, Shipping, and Serving of patient monitor:

Shipment: We are shipping products after receive payment 7 work days

Delivery time: 7 work days

Package: Standard package

After sale services: We provide one year's warranty. We guarantee help customers to settle all of quality problem with free during warranty.

Related Products:



patient monitor YJ-F15

patient monitor YJ-F5



patient monitor YJ-FM2.5



Our Services

1. OEM/ODM service : Forever Medical offers OEM/ODM service for distributor to establish long term business relationship.
2. Warranty : Forever Medical offers 18 months free after-sale service for customers and be responsible for life-time repair.
3. Technical support: Forever Medical offers professional engineers for customers on the technical support, and response within 12 hours

[Your satisfaction is our duty! It will be our great honor to serve you. We do hope we can offer you more help in medical field and build the long term business relationship with you.](#)